

### Method Information

#### Method Comments

Method Modified	User kiran_125274
Method Locked	No
Method	Id 5021
Old Id	
Method Version	3
Method Edit User	kiran_125274
Source S/W Info	Empower 3 Software Build 3471 SPs Installed: Service Release 4 DB ID: 3125586228

### W2690/5 instrument Setup

Type	W2690/5
Instrument Status	On
Channel Name	2690/5
Description	
Use Channel Monitor	Off
Monitor Parameter	System Pressure
Stroke Volume	130uL (flow rates <= 10.000 mL/min)
Chart Out	%A
Syringe Draw Rate	Normal
Depth Of Needle	0.0
Degas Mode	On
Pump Mode	Gradient
Flow	1500
PA	52.0
%B	48.0
%C	0.0
%D	0.0
High Limit	5000.0
Low Limit	0.0
Enable Sample Temp	True
Enable Column Temp	True

Bubble Detect	True
Pre Column Volume	0.0
Sample Temp Target	5.0
Sample Temp Range	2.0
Sparge A	0.0
Sparge B	0.0
Sparge C	0.0
Sparge D	0.0
Column Temp Target	45.07
Column Temp Range	2.0
Flow Ramp	2.00
Column choice	No Change